

NITROGEN GENERATORS



HGT-HN01

- Widely used to supply nitrogen to devices such as GC, GC-MS, GCxGC-MS, DART ion source, etc.
 High-purity nitrogen is suitable for laboratories, with purity comparable to cylinder gas
 The strictly selected molecular sieves can effectively remove oxygen, water, hydrocarbons and other impurities
- 304 stainless steel gas storage tank and nitrogen generation tower

SPECIFICATION

or Edit Iowillow		
Code	HGT-HN01	HGT-HN02
Output gas flow rate	1L/min	2L/min
Output gas pressure	80psi	
Hydrogen purity	99.999%	
Noise	≤ 55dB	
Output interface	1/8 (can be customized)	
Power supply	AC220V, 50Hz, 1.2kW	
Dimensions (L×W×H)	430×620×1080mm	
Net weight	110kg	

STANDARD DELIVERY

Main unit	1pc
Tools	1set

OPTIONAL ACCESSORY

Hydrocarbon Eliminator HGT-HC210



hydrocarbon eliminator (optional)